INSTRUCTIONS FOR USE

P/N SSLW-0100 Light Wand

CE 0297

EC/REP

RMS UK, Ltd.
28 Trinity Road
Nailsea, Somerset BS48 4NU
United Kingdom
TEL: 01275 858891

Sunoptic Technologies®
6018 Bowdendale Avenue
Jacksonville, FL 32216 USA

Customer Service: 904 737-7611
FAX 904 733 4832
Toll Free 877-677-2832

PLEASE READ BEFORE USE

Failure to follow these instructions could render the device unusable and void any warranty.

DEVICE DESCRIPTION/DEVICE CLASSIFICATION

SSLW-0100 Light Wand – Class II US FDA, Class IIa EU

SUNOPTIC SURGICAL® fiberoptic light guides are designed to deliver maximum light when coupled to a medical grade fiberoptic light source. SUNOPTIC SURGICAL® fiberoptic light guides can be used with quartz halogen, metal halide, or xenon light sources.

INTENDED USE

This device is designed to illuminate a surgical site by relaying light from a fiberoptic light source on to the desired site. SUNOPTIC SURGICAL® light guides are medical grade, high transmittance, peak efficiency.

CONTRAINDICATIONS

This device has no contraindications that of which we have knowledge.

WARNINGS & PRECAUTIONS

The user of this product should be thoroughly familiar and trained in use and care of the product. Use caution to treat SUNOPTIC SURGICAL® light guides as you would any fine optical device.

SUNOPTIC SURGICAL® light guides are provided non-sterile and must be sterilized before use. See instructions for cleaning and sterilization

Always inspect the light guide for any evidence of damage prior to use. Pay particular attention to optical surfaces looking for scratches or dings

Use caution to treat SUNOPTIC SURGICAL® light guides as you would any fine optical device

For use with 300 watt xenon light source be sure light source is filtering at least 90% of Infrared Radiation to prevent high heat

The user of this product should be thoroughly familiar and trained in use and care of the product

FIRE HAZARD: Never drape or cover the end of any fiberoptic light guide with anything flammable. Never place the unconnected distal end of a light guide on a drape or other flammable material.

Do not use more than two light guides coupled to independent light sources at once.

Do not bend light guide cable to less than a ¾ inch radius. High intensity light may cause tissue burns. User must position light guide at sufficient distance from tissue to avoid burns.

Do not use automatic cleaning machines as they may cause fiber breakage.

ADVERSE EVENTS

We have no knowledge of any adverse events with this product.

CLEANING AND STERILIZATION FOR FIBEROPTIC LIGHT GUIDES

The cleaning and sterilization methods for the devices were validated and information on them is provided under consideration of ISO 17664.

Fiberoptic light guides are high quality optical devices. They require similar care to that taken for any precision optical component.

During cleaning, care must be exercised handling the light guides to avoid breaking fiber. Do not forcefully bend, kink, or handle roughly.

CLEANING

Cleaning solutions must be prepared according to the manufacturer's instructions and compatible with the material components (silicone, stainless steel, and glass) of the light guide.	
Step 1:	Rinse the light guides in running water to remove all loose or soluble contamination.
Step 2:	With fresh, non-protein fixing, cleaning solution, hand wash the Light guides.
Step 3:	Using a soft, non-metallic, brush and lint free towels remove all soil from the light guide. Special attention should be paid to the junctions between flexible and solid parts to remove any trapped soil.
Step 4:	Once clean, thoroughly rinse in running water to remove all traces of the cleaning solution
Step 5:	Rinse in distilled water to remove any dissolved minerals in the cleaning water.
Step 6:	Clean the fiber ends by gently rubbing with a cotton swab saturated in alcohol. Immediately dry using compressed air.

STERILIZATION - STEAM AUTOCLAVE

Sterilization can be accomplished in an AAMI ST8 compliant steam sterilizer with the following criteria:

Preparation: SPSmedical Self-Seal pouch SSP387 with SPSmedical STEAM*Plus* Integrator SSI-100

METHOD	CYCLE	MIN. TEMP.	EXPOSURE TIME	DRY TIME
Steam (in pouch)	Prevacuum	Prevacuum 270 deg. F (132 deg. C)	4 minutes	40 minutes

The above steam sterilization method was validated for these devices. If other equipment, methods or parameters are used, the user is responsible for validation.

LIMITED WARRANTY

Your fiberoptic light guide has a one (1) year warranty from the date of shipment for defects in material/workmanship, **except** for broken fiber. Should your product prove to have such defects within one (1) year of shipment, SUNOPTIC TECHNOLOGIES® will repair or replace the product or component part without charge. Warranty **does not** cover equipment subject to misuse, accidental damage, and normal wear and tear. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state. Should your fiberoptic light guide need servicing under this warranty, please contact your distributor or your customer support specialist for return authorization documentation. You should carefully pack the product in a sturdy carton, including a note describing the defects, your name, your company name, telephone number and a return address.

POST WARRANTY REPAIRS

Please contact your distributor or your customer support specialist for return authorization documentation.

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