HEADLIGHT CABLES FEATURES, CARE AND USE

Sunoptics Surgical[®] headlight cables manufactured by **Sunoptic Technologies[®]** have several important features to increase performance and durability.

FEATURES	DETAILS
The fibers are mounted on the proximal end by a process called “FUSING.”	The fusing serves to fix all fibers together without the use of epoxy. This <i>eliminates</i> the problem of rapidly deteriorating performance or failure due to burning epoxy.
	The process also compacts the fiber and allows a 15% more efficient conductor of light
The cables have a lens mounted in the proximal “WOLF” connector.	The lens is mounted inside the cable to provide a very accurate light focus. This ensures that Sunoptics Surgical[®] cables can transmit the same amount of light through a 3.5mm cable that is normally transmitted through a 5mm cable. Thus, the headlight weight is reduced by 50% and more comfortable for the surgeon.
Cable has a mono-coil embedded in the silicone sleeve	Bifurcated cables have the mono-coil embedded in the sleeve from the light source connector to the junction block. This is the area most abused through handling. The branches are left without mono-coil to ensure light weight and comfort.
Cables have “Ultra Grip” handle	The “Ultra Grip” handle is a design feature to ensure that the user removes the cable from the light source using the “Ultra Grip” handle instead of removing by pulling on the sleeve. The removal by pulling on the sleeve will cause the sleeve to detach from the handle and reduce the cable life.
	The “Ultra Grip” handle fixes the angle the fiber bends in a rigid environment, reducing repeated bending when attaching and detaching the cable to the light source. It also ensures the fiber cable will point toward the floor without introducing an additional stress point.
CARE AND USE	
Fiberoptic cables are manufactured with at least 5,000 ultra-small fibers, manufactured from glass. Therefore, even though we add many protective features, the heart of the cable is still manufactured using glass.	When removing and storing the cable ensure that the diameter of the cable coil is at least 12 inches (30cm’s).
Cable fusing increases heat durability and performance but the process does create a short length of solid glass.	This can be cracked like any piece of glass and needs to be handled with care, for example, a blow such as falling on a hard floor can crack the glass.
The cable lens mounted in the light source connector is a very important performance enhancement to the cable.	It is also a piece of glass and may be cracked by a sudden sharp blow such as hitting the hard floor.
Cleaning	The cable is NOT intended to be autoclaved or sterilized using liquid sterilizing chemicals. Cleaning can be accomplished using whatever cleaning agents are used in the hospital.

SUNOPTIC TECHNOLOGIES[®]
6018 BOWDENDALE AVENUE •JACKSONVILLE, FLORIDA 32216 USA
PHONE : 904-737-7611 •FAX : 904-733-4832

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