LE-3000

LED Light Engine

Installation and Use Manual

Sunoptic Technologies®
6018 Bowdendale Avenue
Jacksonville, FL 32216 USA

Customer Service: 904 737 7611
Toll Free 877-677-2832

RMS UK, Ltd.
28 Trinity Road
Nailsea, Somerset BS48 4NU
United Kingdom
TEL: 01275 858891

TABLE OF CONTENTS

1. INTRODUCTION
2. WARNINGS
3. INSTALLATION
4. CLEANING
5. DISINFECTION
6. SERVICING & REPAIR
7. END OF PRODUCT LIFE
8. CHART OF MEDICAL DEVICE SYMBOLS
9. APPENDIX A

1. INTRODUCTION

Congratulations on the purchase of your new LE-3000 Light Engine!

This high efficiency light engine uses state-of-the-art superior illumination technology and offers a variety of features such as:

- 7000^o K daylight brightness for good color definition
- Quiet operation
- Long life, 30,000 hrs.
- Compact and light weight
- Turret which adapts to various types of light guides

In short, you have chosen the best and we would like to make sure you receive the optimal results with your new LE-3000 Light Engine by using it correctly.

This Installation and Use Manual will help you to install the device and optimally integrate it with other components of your system.

2. WARNINGS/ CAUTION

 **Caution** To prevent fire or electric shock, do not open or expose the LE-3000 Light Engine unit to rain or moisture. Refer all servicing to qualified personnel only.

 **Caution** Not suitable for use in presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide.

 **Caution** To prevent any potential electro-magnetic interference, do not use any kind of cellular phone near the light engine.

 **Caution** This product should be used only with type BF endoscopic instruments which have been certified according to IEC 60601-1 for medical equipment and IEC 60101-2-18 for endoscopic equipment.

 **WARNING**
Do not modify this equipment without authorization of the manufacturer.

 **Caution** All devices connecting to the LE-3000 Light Engine may be classified as medical equipment. When additional information processing equipment is connected to the LE-3000 Light Engine, the operator must determine that all equipment complies with the appropriate end-product standards (such as IEC 60950 or IEC 60065 and the Standard for Medical System, IEC 60601-1-1).

 **WARNING** The LE 3000 is a highly concentrated light engine (luminous power per area) and this high energy density is retained through connected lightguides and instruments. The output of a connected instrument left in close proximity to tissue or flammable materials presents a risk of injury or fire. Qualified personnel must determine a safe working distance and intensity setting for each application. The output should never be left on unattended. Turn the light engine off or place it in standby if it will not be required for a period of time.

 **Caution** The LE-3000 Light Engine when installed into light source can cause permanent eye damage if viewed directly with unprotected eye. To reduce the chance of eye damage, set the intensity control **always** to the minimum level and plug the fiber optic cable into the unit before turning **on** the power.

EQUIPMENT CONNECTION

 **Caution** The fiber optic cable should be a **NON-CONDUCTIVE CABLE** unless otherwise determined by the end-user integrating this device into an end product. As such, attached equipment should not have conductive shielding or any other conductive connection between the patient and equipment. Such connection will impair safety of the equipment. It must be rinsed free of soaking/disinfectant solution and dried before plugging into the LE-3000 Light Engine receptacle. Ensure the optical surface is clean before engaging into the light source.

 **Caution** The proximal end tip of installed lightguides and accessories can become hot. Allow end tips to cool before removing.

3. INSTALLATION

NOTE Avoid places where the light engine may be splashed with liquid.

Absolutely DO NOT use in any environment with explosive or flammable gases.

CONNECTING THE FIBER OPTIC CABLE

Connect the light cable to the instrument and then plug the light cable end fitting into the turret port on the front panel.

NOTE The turret has a “**cable in place**” sensing feature that will not allow the LED light to turn on without a fiber optic cable inserted into the turret port.

4. CLEANING

CLEANING THE UNIT

The exposed turret can be cleaned with any cleaning agent, which is used for external cleaning of electric equipment and according to instructions given by the manufacturer of the cleaning solution.

Do not allow moisture or liquids to come in direct contact with the unit except for those used in cleaning the exposed turret.

Do not use cleaning agents on the exposed turret that are not permitted for use with plastics, i.e., ammonia, acetone, salty acids (HCl), etc.

CLEANING THE OPTICS

Exercising proper care for the glass optic in your LE-3000 Light Engine can help maintain light quality and extend usable lifetime. Using the proper cleaning products and methods are as important as cleaning the component itself.

Cleaning Solutions for the Optics

There are a variety of cleaning solutions you can use such as Reagent-Grade Isopropyl Alcohol, Reagent-Grade Acetone and De-Ionized Water to ensure a long product lifetime. A good solvent mix to use is a mix of 60% Acetone* and 40% Methanol. (Acetone alone dries too quickly to dissolve all of the debris. The Methanol slows the evaporation time, and also dissolves debris that Acetone alone would not clean.)

NOTE Acetone should never be used on plastic optics or components as it will damage them.

Cleaning Procedure

Step #1

First use compressed air (aerosol can) to remove surface dust without directly contacting any coating the optic may have. Wiping a dusty optic with cotton swabs or lens paper is like cleaning it with sandpaper.

Step #2

Soak with solvent a Cotton-Tipped Swabs or Lens Tissue then use a circular motion starting from the edge of the optic moving towards the center to wipe away any dust or dirt. Repeat as necessary; each time use a new cotton swab or lens tissue.

5. DISINFECTION

DISINFECTING THE EXPOSED TURRET ONLY

Use any disinfectant agents which are commonly applied while disinfecting surfaces of electric medical equipment. Such disinfectant agents are usually in the form of sprays or damp cloths.

Follow the instructions given by the manufacturer of the disinfectant solution.

6. SERVICING AND REPAIR

Defective items of equipment are to be serviced and repaired exclusively by persons authorized by the manufacturer. All repair work shall employ original manufacturer's parts only.

LIMITED WARRANTY

Your LE-3000 Light Engine carries a 2-year limited warranty from the date of shipment on workmanship and all defects of material. Should your product prove to have such defects within two years from date of shipment from **Sunoptic Technologies®**, we will repair or replace the product or component part without charge. Should your LE-3000 Light Engine need servicing under this warranty, please contact **Sunoptic Technologies®** for return authorization documentation. You should carefully pack the unit in a sturdy carton and ship it to the factory. Please include a note describing the defects, your name, telephone number and a return address. Warranty does not cover equipment subject to misuse, accidental damage, normal wear and tear or if transferred to a new owner without authorization from **Sunoptic Technologies®**. This warranty gives you specific legal rights and you may also have other rights that vary from state to state.

POST WARRANTY REPAIRS: You may return your product(s) for repair, shipping prepaid to the factory. Your product will be inspected and an estimate of repair charges will be submitted to you for approval. Payment must be received before repairs are completed.

- In the US: 877-677-2832 (toll free)
- FAX number: 904 733-4832
- Customer Service 904-737-7611

7. END OF PRODUCT LIFE

In accordance with the European Waste from Electrical and Electronic Equipment (WEEE) directive, we encourage our customers to recycle this product whenever possible. Disposal of this unit must be performed in accordance with the applicable local environmental regulations.

In the US a list of recyclers in your area can be found at: <http://www.eiae.org/>.

Please contact customer service to issue a return authorization to return product to manufacturer at the end of product life.

8. CHART OF MEDICAL DEVICE SYMBOLS USED

	Manufacturer
	Date of manufacture (YYYYMM or YYYY)
	Symbol for "Authorized Representative in the European Community."
	Caution, consult accompanying documents
	CE mark


9.0 APPENDIX

LE-3000 LED Light Engine Specifications

(due to continuous product improvement, subject to change without notice)

Optical Performance

Luminous flux through a 5mm aperture into 0.66NA	≥1400 Lumens or greater
Luminous flux through a 5mm aperture into 1.0NA	1900 Lumens nominal
Color temperature	7000K nominal
Color rendering index	78 nominal
Lamp life	30000 hours

Lamp Life; Dependent on operating junction temperature; controlled by end-user

Electrical specifications

- Input power: 12VDC, 60 – 90W continuous
- Power to LED: 50 – 75W continuous
- LED drive current: Fixed 18A_{RMS} LED current or user configurable 4 – 18A_{RMS}
- Electrical efficiency: typical $\geq 83\%$
- Maximum operating frequency: 500kHz

Absolute Maximum Ratings	
Local ambient operating temperature airflow	-10 to 40°C with ≥ 20 CFM
Storage temperature	-40°C to 85°C
LED junction temperature	150°C
PCB temperature	85°C
Input power voltage	13.2 VDC
Maximum accessory current draw	500mA
Voltage on VDD with respect to DGND	2.7V to 5.75V
Voltage on ADJ pin with respect to PS_POS	-0.3V to 7.0V
Voltage on EN, DIM with respect to DGND VDD+0.5V	-0.5V to
Voltage on END_TIP with respect to DGND	-5V to 30V


Input Power Connector (J1)					
Header: 0.156" Friction lock, 4 position					
Molex PN: 09-65-2048					
Pin	Description	Conditions	Min	Typ	Max
1	POSITIVE	Pins 1&2 connected	10.8V	12V	13.2V
2	POSITIVE		10.8V	12V	13.2V
3	NEGATIVE	Pins 3&4 connected	-	-	-
4	NEGATIVE		-	-	-

Control Connector (J2)					
Header: 0.1" Friction lock, 12 position					
Molex PN: 22-23-2121					
Pin	Description	Conditions	Min	Typ	Max
1	ADJ – LED current reference voltage	With respect to CGND	0.3V	-	1.5V
2	EN – Input, LED drive enable	Logic level limits	≤0.8V (LOW)	-	≥2.0V (HIGH)
3	DIM – Input, PWM dimming signal	Logic level limits	≤0.8V (LOW)	-	≥2.0V (HIGH)
4	LED_TEMP – Output, LED temperature sensor	Thermistor referenced to DGND	-	10kΩ @25°C	-
5	PCB_TEMP – Output, PCB temperature sensor	Thermistor referenced to DGND	-	10kΩ @25°C	-
6	END_TIP – Output, Fiber optic cable end tip sensor	Open collector to DGND. Recommend 10k pull up.	-	-	5mA
7	FAN_CTRL – Input, Fan enable	Optoisolator LED input. Requires limiting resistor.	7mA	10mA	13mA
8	FAN_TACH – Output, fan tachometer	Open collector to DGND. Recommend 10k pull up	-	-	5mA
9	+5V – Output, 5V accessory power	With respect to CGND Maximum draw 500mA	4.95V	-	5.15V
10	VDD – Input, user positive reference	With respect to DGND. Short JP2 to connect to +5V	2.7V	-	5.75V
11	CGND – Output, ADJ and +5V ground reference	Connected to input power supply positive terminal.	-	-	-
12	DGND – Input, user ground reference	Short JP3 to connect to CGND	-	-	-


External Potentiometer Connector (J5)					
Header: 0.100" Friction lock, 3 position					
Molex PN: 22-23-2031					
Pin	Description	Conditions	Min	Typ	Max
1	Main Terminal	High side	-	10kΩ	-
2	Wiper		-	-	-
3	Main Terminal	Low side	-	-	-

Mechanical Specifications


(due to continuous product improvement, subject to change without notice)

3.60 HEIGHT X 4.83 WIDTH X 6.95 DEPTH
 2.7 LBS WITH TURRET
 2.05 LBS WITHOUT TURRET

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