# LED ROVER IV

## Instructions for Use

**DIST**

mytaMed  
274 Lancaster Ave., Suite 208  
Malvern, PA 19355  
Tel: 877 444 6982  
www.mytamed.com

Sunoptic Technologies®  
6018 Bowdendale Avenue  
Jacksonville, FL 32216  
Customer Service: 904 737 7611  
Toll Free 877 677 2832

**EC|REP**

RMS UK, Ltd.  
28 Trinity Road  
Nailsea, Somerset BS48 4NU  
United Kingdom  
TEL: 01275 858891

## TABLE OF CONTENTS

1. INTRODUCTION
2. INTENDED USE/INDICATIONS FOR USE
3. WARNINGS
4. SPECIFICATIONS
5. OPERATION
6. BATTERY INFORMATION
7. CLEANING/STERILIZATION
8. WARRANTY
9. CHART OF MEDICAL DEVICE SYMBOLS

## 1. INTRODUCTION

The LED Rover IV is a Light Emitting Diode (LED) and battery unit designed to couple onto the light post of commercially available endoscopes having the common arrangement of an ACMI post with nesting screw-on adaptors for Wolf and Storz / Olympus. The LED Rover IV is made to take the place of conventional remote light source and fiberoptic cable.

To connect it to the endoscope remove both adaptors if they are present. If your mytaMed Rover IV has an ACMI port, it couples to the endoscope in the same fashion as a conventional light cable by press fitting onto the ACMI post. If your mytaMed Rover IV has a threaded port, it screws onto the threads for the Wolf adaptor.

The LED Rover IV is powered by a single "18650" size lithium-ion 3.7V rechargeable battery.

The battery powers an LED via a current converter that supplies a constant current to the LED. The light is directed into the light post of the scope. The light is then carried by fiber optics internal to the scope to illuminate the target.

## 2. INTENDED USE / INDICATIONS FOR USE

The LED Rover IV is used to illuminate the site of surgery during minimally invasive surgical procedures in arthroscopy (orthopedic surgery), laparoscopy (general and gynecological surgery) and in Endoscopy (general, gastroenterological and ENT surgery). The light is typically transmitted from the light source through a scope.

### 3. WARNINGS

#### 3.1 WARNINGS

##### CAUTION

Federal law restricts this device to sale by or on the order of a licensed practitioner.

##### CAUTION

Battery powered equipment is considered unsuitable for use in the presence of flammable anesthetics.

##### CAUTION

The LED Rover IV is provided **non-sterile**. If used for sterile procedures, the unit must be sterilized or draped with commercially available surgical drapes.

##### CAUTION

The LED Rover IV will become warm to the touch but even prolonged contact during use will not cause injury.

##### CAUTION

Avoid looking directly at the LED inside the port or pointing it directly at anyone's eyes. The LED is painfully bright, can momentarily blind a person and can cause lingering ghost images.

##### CAUTION

Handle damaged or leaking batteries with care. Should you come in contact with the electrolyte, wash the exposed area with soap and water. Should the eye come in contact with the electrolyte, flush the eye with water for 15 minutes and seek medical attention. Do not incinerate discarded batteries. Follow local EPA requirements for proper disposal.

##### CAUTION

Always follow the instructions for use provided with the rechargeable batteries and charger system which are provisioned for use with your ROVER IV system.

##### CAUTION

Only operate the Rover IV with batteries and battery charger approved by the manufacturer.

#### 3.2 CONTRAINDICATIONS

***The device has no known contraindications.***

## 4. SPECIFICATIONS

Item	Specification
Lamp Type	LED
Color temperature	5000° K (nominal)
Lamp life	50,000 hrs. (nominal)
Battery	Lithium Ion 18650 3.7V
Current converter output	611 milliamps (nominal)
Endurance time	2.5 hours at full output, steady decline for several hours afterward
Model	LED Rover IV
Dimensions	5.0 inches long X 1.44 inches diameter at LED housing
Weight	5.5 oz.
Electrical Protection	Class I, Type BF
Medical Classification	FDA Class II Europe Class I

## 5. OPERATION

The LED Rover IV is operated by a pushbutton switch inside a protective silicone rubber boot at the back end of the handle. It features a momentary ON function with partial depression of the switch, and clicks to stay ON with full depression. Depress the switch fully again to turn it off.

## 6. BATTERY INFORMATION

### 6.1 RECHARGEABLE BATTERY

The rechargeable batteries provided are high capacity 18650 3.7 volt Lithium-ion battery. The NITECORE Intellicarger i2 battery charger provided is compatible only with the 3.7 volt high capacity rechargeable Lithium-ion batteries that are supplied. The 3.7 volt high capacity batteries provided are compatible *only* with chargers specifically designed for them. Do not put attempt to charge the batteries with another type of charger.

A fully-discharged battery may require approximately eight hours to reach a full charge. The three indicator lights on the charging cradle will be blank when there is no battery in the cradle. The charger will flash one LED while determining what configuration of battery is in the charger and this LED will remain solid once the configuration is determined. The second and third LEDs will flash while it is charging, and there will be three solid lights when the battery is fully charged. Batteries may be "topped off" to full charge even if they are not fully discharged. Batteries may be left in the charging cradle when charged. Avoid directly short circuiting a battery. The power converter operates on 100~240 VAC, 50~60Hz @ 0.25A.

With fully charged batteries the light source is approximately 2.5 hours with the strongest brightness. The following 5 hours there will be a slow light waste.

**Consult the charger's individual instructions for further operational details, additional warnings, and safety data.**

## 6.2 BATTERY REPLACEMENT

Unscrew the headpiece that houses the LED. (Do not remove the switch; the battery cannot fit through the switch end.) Remove the discharged battery and insert a fully charged battery. Observe correct polarity. The positive (+) end of the battery must point toward the headpiece. The positive end of the battery has a protruding nipple; the negative (-) end is flat.

## 7. CLEANING / STERILIZATION

### 7.1 CLEANING

The LED Rover IV is a high quality optical device. Care for it appropriately.

A soft brush with mild detergent may be used to clean the device. The lens of the LED inside the optic port may be wiped with a swab. Avoid harsh or abrasive material that could scratch or in any way damage the optical surface.

### 7.2 STERILIZATION

The LED Rover IV is provided non-sterile. The device should be cleaned and disinfected or sterilized prior to use, according to your institution's requirements.

#### **CAUTION**

***Do not reprocess the LED Rover IV using steam sterilization or autoclave. Use of these processes will result in damage to the instrument, and void its warranty.***

Process per your institution's validated procedures and process parameters.

Ensure that fresh or adequately charged batteries are installed and that the headpiece and tailpiece are snugly tightened before sterilizing or disinfecting.

**Compatible Sterilization Method: STERRAD™ in STERRAD sterilizer.** Double wrap with sterilization wrap. The full cycle of the Sterrad sterilizer comprises two identical exposure phases (6 minute injection, 2 minute diffusion and 2 minute plasma) preceded by a 10-minute pre-exposure plasma phase.

Follow the instructions provided with the STERRAD Sterilization System for sterilization processing via this equipment and method.

**NOTE: The LED Rover IV has a black anodized exterior finish which is affected by Sterrad sterilization. The black anodized finish gradually bleaches; however, this is irrelevant to the functionality of the unit.**

### 7.3 HIGH LEVEL DISINFECTION (IMMERSION)

The LED Rover IV may be immersed in high level disinfectant solution no deeper than six (6) inches (15cm) for as long as one (1) hour. Process according to the disinfectant manufacturer's instructions or your institution's established protocol. We make no claims other than the Rover's IV immersibility as stated.

## 8. LIMITED WARRANTY

Your LED Rover IV carry a 1-year warranty from the date of shipment on workmanship and all defects of material except batteries. Should your product prove to have such defects within one year from the date of shipment, **mytaMed** will repair or replace the product or component part without charge. Should your fiberoptic product(s) need servicing under this warranty, please contact **mytaMed** for return authorization documentation. You should carefully pack product in a sturdy carton and ship it to the factory. Please include a note describing the defects, your name, telephone number and a return address. Warranty does not cover equipment subject to misuse, accidental damage, and normal wear and tear. This warranty gives you specific legal rights and you may also have other rights that vary from state to state.

- In the U.S call: 877 444 6982 (toll free)
- E-mail: [customerservice@mytamed.com](mailto:customerservice@mytamed.com)

## 9. CHART OF MEDICAL DEVICE SYMBOLS USED

	Manufacturer
	Date of manufacture (YYYYMM or YYYY)
	"Authorized Representative" in the European Community
	Caution, consult accompanying documents
	Consult Instructions for Use
	CE mark

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