



SSLP-050T 50W LED Light Source with Headlight System in Carrying Case



Portable LED Light Source with Headlight System

- 50W medical grade light source with 4-port turret (Acmi, Wolf, Storz, and Olympus)
- Color temperature: 6500° K
- Typical LED life: >50,000 hours
- Accepts 100-240V AC, 50/60 Hz input power for worldwide operation
- Deluxe bifurcated headband with top and rear knob adjustment
- Headlight module with variable spot adjustment 20-80mm
- 5mm bifurcated headlight cable for optimum light output
- Universal clamp mount and bracket for pole mounting LED light source included
- Durable carrying case to safely store and transport complete system
- Additional cutout for small scope (scope not provided)
- CE marked

Sunoptic Technologies®, 6018 Bowdendale Avenue, Jacksonville, FL 32216
Toll Free: 877-677-2832 Fax: 904-733-4832

www.sunoptictech.com

